

microphone instruments



ReSound

rediscover hearing

About ReSound

ReSound provides excellent sound by offering innovative hearing solutions that combine original thinking and design with solid technology — all based on deep audiological insight and a very good understanding of users.

For more information about **ReSound Live™** contact your hearing healthcare professional today.

www.resound.com/Live

ReSound and rediscover hearing are registered trademarks of The GN ReSound Group. © 2009 The GN ReSound Group, all rights reserved.



M101070.GB.09.11. Rev. A

*Source: 2008, trials conducted internally at ReSound

ReSound Live™

Who's smiling now?



You'll be smiling

when you wear **ReSound Live™** remote



ReSound

rediscover hearing

A revolution in your ear

New remote microphone technology (RM) instruments might be just what you need

New technology for people just like you

You want your hearing difficulties solved by a small, discreet instrument hidden away in your ear. But until recently, if you had more severe hearing loss, a different inner ear shape or narrow ear canals, it was difficult to provide enough hearing assistance without increasing the size — and obviousness — of the instrument.

Now ReSound's remote microphone technology (RM) has changed all that. In fact, it's probably the most significant development in in-the-ear instruments in the last 20 years. Thanks to their entirely new design, the new ReSound Live™ family of RM instruments delivers:

- Smaller size
- More power
- Better sound quality

RM instruments are perfect for anyone who needs more power or battery life — or who just wish their current in-the-ear, in-the-canal or completely-in-the-canal instruments were a little less obvious.

Much more than you can see

The first real innovation in custom instruments in 20 years



Conventional in-the-ear instruments vs. significantly smaller Remote Microphone instruments

What a difference! Remote microphone technology takes advantage of the ear's anatomy to place the microphone in the area known as the 'concha'. This reduces instrument size and lifts performance, providing more power and battery life in a less visible and more secure fitting than with any other product worn in the ear.



More than just good looks

ReSound Live RM instruments do more than just hide in your ear. They apply the best of today's technology to help you hear better in all kinds of situations. But that's not all. The unique positioning of the microphone makes the most of the ear's naturally sophisticated acoustics to deliver stronger performance in speech understanding, awareness, and wind noise protection. So listening becomes much more comfortable and natural.

Natural directivity

Your ear is shaped to help you prioritise voices and pinpoint where sounds are coming from. RM instruments use this natural design to help you participate in conversations and maintain awareness of your surroundings.

Natural, undisturbed sound quality

Every ReSound Live™ instrument features a revolutionary sound processor that delivers richer, more detailed and more true-to-life sound quality. Additionally, WhistleControl™ cancels whistling without sacrificing sound quality or audibility.

Natural wind noise protection

Placing the microphone within the contours of the ear assists to minimise wind noise — so much so that in a recent test of RM instruments, 97% stated they weren't bothered by wind noise*.



Feature Highlights

- Smaller size even for severe hearing loss
- Natural directivity for better speech understanding
- Natural localisation of sounds around you
- Natural wind noise protection
- Dual Stabilizer® II DFS with WhistleControl™ feedback cancellation
- NoiseTracker™ II noise reduction for greater comfort
- 17-band Warp™ compression for clear, undistorted sound processing
- Up to 176 hours of battery life
- PhoneNow™ (optional) for improved phone usage

