

*Quick Reference:*

**YOUR CUSTOM  
COMPLETELY-IN-THE-CANAL  
(CIC) REMOTE MICROPHONE  
HEARING INSTRUMENTS**

*Model information:*

Your hearing instrument model designation  
is \_\_\_\_\_.

*Battery Information:*

Your hearing instrument uses size  
10A battery.

## **CUSTOM COMPLETELY-IN-THE-CANAL (CIC) REMOTE MICROPHONE HEARING INSTRUMENTS**

*This booklet contains information regarding models: LV10-M, LV10-MP, be910-M, be910-MP, be710-M, be710-MP.*

1. Body
2. Battery compartment and on/off switch
3. Sound outlet
4. Wax filter
5. Microphone tubing
6. Microphone
7. Microphone sound inlets
8. Air vent
9. Manufacturer, model and serial number (imprinted on side of body)

