

ReSound



Quick Reference:

**YOUR INVISIBLE OPEN
TECHNOLOGY (IOT)
REMOTE MICROPHONE
HEARING INSTRUMENTS**

Model information:

Your hearing instrument model designation
is _____.

Battery Information:


Your hearing instrument uses size
10A battery.

MK602124 Rev.A

10 IOT

INVISIBLE OPEN TECHNOLOGY (IOT) REMOTE MICROPHONE HEARING INSTRUMENTS

*This booklet contains information
regarding models: be900, be700.*



1. Body
2. Battery compartment and on/off switch
3. Sound outlet
4. Microphone tubing
5. Microphone
6. Microphone sound inlets
7. Serial number
8. Manufacturer and model
9. Wax filter replacement kit

