

People with hearing loss are at the heart of everything we do. Our award-winning ReSound Smart Hearing aids help people feel more involved, connected and in control. In an ever-smarter world, we think big and challenge the norm so that we can transform lives through the power of sound. Available exclusively through audiologists and hearing care professionals.

ReSound® is part of The GN Group - pioneering great sound from world-leading ReSound hearing aids to Jabra office headsets and sports headphones. Founded in 1869, employing over 5,000 people, and listed on NSDAQ OMX Copenhagen, GN makes life sound better.



### [resound.com/wirelessconnectivity](https://resound.com/wirelessconnectivity)

© 2016 GN ReSound A/S. All rights reserved. ReSound is a trademark of GN ReSound A/S. Apple, the Apple logo, iPhone, iPad, iPod, iPod touch, FaceTime and Apple Watch are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc. Android and Google Play are trademarks of Google Inc. The Android robot is reproduced or modified from work created and shared by Google and used according to terms described in the Creative Commons 3.0 Attribution License. Dolby and the double-D symbol are registered trademarks of Dolby Laboratories. Bluetooth is a trademark of Bluetooth SIG, Inc.

#### Manufacturer according to FDA

#### RESOUND NORTH AMERICA

8001 E. Bloomington Freeway  
Bloomington, MN 55420  
1.888.735.4327  
resound.com

#### RESOUND GOVERNMENT SERVICES

8001 E. Bloomington Freeway  
Bloomington, MN 55420  
1.800.392.9932  
resound.com/veterans

#### Manufacturer according to Health Canada

#### RESOUND CANADA

303 Supertest Road  
Toronto, Ontario M3J 2M4  
1.888.737.6863  
resound.com



# Expand your hearing experience

SMART HEARING ECOSYSTEM OF WIRELESS CONNECTIVITY



MK604656 Rev C 2016.11



GN Making Life Sound Better

# Connect, interact and engage with your world

With our Smart Hearing™ ecosystem of wireless connectivity you can stream, control and personalize the sound to your ReSound hearing aids without wearing an intermediary device around your neck.



**MADE FOR IPHONE™**  
Make your hearing aids work like wireless stereo headphones.



**WIRELESS ACCESSORIES**  
Enjoy opportunities to hear and do more in challenging listening situations.



**APPS**  
Control and personalize your hearing experience from your iPhone or Android™ smartphone.



**RESOUND MICRO MIC**  
Focus on the sounds you want to hear, even in a noisy place, and extend your hearing range to over 80 feet.



**RESOUND MULTI MIC**  
Get all the benefits of the ReSound Micro Mic, and more. When placed on a table, the ReSound Multi Mic automatically turns into a table microphone. Plus, it connects with loop and FM systems\*, and uses a line-in for streaming audio to your hearing aids.



**RESOUND PHONE CLIP+**  
Make crystal-clear phone calls, mute background noise while you talk, or stream music, podcasts or any audio from your phone.



**RESOUND REMOTE CONTROL 2**  
Adjust the volume or mute your hearing aids, switch programs, and see all your settings at a glance.



**RESOUND TV STREAMER 2**  
Stream stereo sound from your TV, computer or music system.



**RESOUND SMART™ APP AND RESOUND CONTROL™ APP**  
Personalize and control your hearing aid settings from your iPhone or Android smartphone.

\* FM receiver required



# Make your hearing aids work like wireless headphones



Made for  
iPod iPhone iPad

Make a phone call, listen to a movie or get turn-by-turn directions all streamed from your iPhone®, iPad® or iPod® touch directly to your hearing aids. With no neck loop or intermediary device, Made for iPhone connectivity is the easiest, most discreet way to make your hearing aids work like wireless stereo headphones.

If you don't have an iPhone, you can still stream stereo sound for calls and music from your Android™ smartphone or any other Bluetooth-enabled phone. Simply pair your phone with the ReSound Phone Clip+.



# Extend your hearing

The latest additions to our range of wireless accessories, the ReSound Micro Mic and ReSound Multi Mic, enable you to connect to the world's smartest hearing aids\*, and feel confident that you won't miss a beat in conversation.



## RESOUND MICRO MIC

Extremely easy to use, the ReSound Micro Mic pairs to ReSound Smart Hearing aids\* in seconds. Clip it on to the clothing of the person you want to talk to and enjoy the conversation, even in background noise – and you can still hear them even if they are over 80 feet away.



## RESOUND MULTI MIC

Works like the ReSound Micro Mic, with the same hearing range of over 80 feet, but also doubles as a table microphone, so you can hear everyone sitting around it. It's compatible with loop and FM systems\*\*, and even has a line-in so you can stream virtually any sound directly to your hearing aids.

\*Compatible with all ReSound wireless hearing aids

\*\*FM receiver required



# Connect with ReSound wireless accessories

## STREAM STEREO SOUND TO YOUR HEARING AIDS

Don't have an iPhone? You can still stream stereo sound to your ReSound hearing aids. Simply pair any Bluetooth-enabled phone with the ReSound Phone Clip+ and use your hearing aids as headphones for phone calls, or for streaming stereo sound from videos, playlists or podcasts.

The Phone Clip+ also works as a basic remote control for your hearing aids. You can adjust the volume and mute/unmute the sound you stream to your hearing aids.



ReSound TV Streamer 2

## ENJOY FAMILY MOVIE NIGHT

With the ReSound TV Streamer 2 plugged into your TV, stereo or computer, you can stream stereo sound directly into your ReSound hearing aids. You can listen at your volume while the people around you listen at theirs.

And, because you can adjust the volume of your surroundings and the TV independently, you can still chat with family and friends while you watch TV together.



Wireless connectivity with the ReSound Phone Clip+.



ReSound Remote Control 2

## CONTROL AT YOUR FINGERTIPS

The ReSound Remote Control 2 gives you simple control and an overview of the volume and program settings of your hearing aids. You can switch to one of your ReSound streaming accessories, balance its volume with your hearing aids and see all your settings on the backlit display.

# Personalize your hearing experience



With the award-winning ReSound Smart™ app, you can fine-tune your ReSound Smart Hearing aids to your unique taste - all without having to touch your hearing aids.

## MAKE ADJUSTMENTS FOR LISTENING COMFORT

Fine-tune hearing aid settings with the Sound Enhancer for more listening comfort in a noisy place or outside on a windy day. You can adjust volume, treble and bass and create favorite settings, geotag them, and have them activate as you enter a new location. And you can control the volume of your ReSound streaming accessories.

## CONTROL THEM FROM YOUR WRIST

With the ReSound Smart app for Apple Watch, you can control your hearing aids straight from your wrist - wherever you are.



## DISCREET CONTROL OF YOUR HEARING AIDS

Adjusting volume or switching programs is as easy as a tap and a swipe. Simply pair your smartphone with the ReSound Phone Clip+ and use the free ReSound Control app for discreet control over your hearing aids. You can also access and balance the volume of your ReSound streaming accessories such as the ReSound Micro Mic.

