



ReSound GN

ReSound LiNX<sup>2</sup>

Hearing has never  
been so effortless



GN Making Life Sound Better



## The latest in Smart Hearing™

Put on your Resound LiNX<sup>2</sup>™ hearing aids and discover the benefits of Smart Hearing that adapts to you and your lifestyle.

Enjoy top-rated sound quality and listening comfort from sunrise to sunset. Get a clear sense of the space and sounds around you with Spatial Sense™.

When you get to work or meet up with friends, engage in conversations knowing you'll follow every word effortlessly, even when noise levels get high – with the help of Binaural Directionality™ II.

Take a phone call as you pull on your jacket and walk across the windy parking lot. As you drive home, stream your latest playlist, audiobook or podcast from your iPhone® without wires or something dangling from your neck.

At the end of the day, you'll still have energy to spare for your loved ones. Because all that hearing, which used to wear you out, now seems so natural and effortless.

Get a vivid  
sense of space  
and speech

**Binaural Directionality II with Spatial Sense** is a combination of formulas unique to Surround Sound by ReSound™ – our sound processing strategy. Your ReSound LiNX<sup>2</sup> hearing aids continuously exchange data to deliver vivid awareness of your sound environment, top-rated sound quality and clear speech understanding, no matter how complex the listening situation.

**Spatial Sense** emulates the way your outer ear channels sound to your inner ear, supporting your ability to detect where sounds are coming from so you experience a more natural sense of your surroundings – with the highest possible sound quality.

**Binaural Directionality II** helps you zoom in on the sounds you're most interested in, while maintaining a sense of what is going on around you. This supports your natural ability to focus on speech, even in the noisiest situations.

This is the essence of ReSound LiNX<sup>2</sup>. No compromise in sound quality or listening comfort. Just effortless hearing as nature intended.



AUDIOLOGY

## Experience Surround Sound by ReSound

Inspired by nature, Surround Sound by ReSound is our unique system for processing sound. Using advanced technologies that emulate the human ear, we take full advantage of your brain's natural ability to recognize sounds and pinpoint where they are coming from.

In a recent independent benchmark test of the latest hearing aids from leading manufacturers, experienced hearing aid users gave Surround Sound by ReSound a top rating across a variety of different sound environments.\*



\* 2013 Benchmark test, DELTA SenseLab.

# Personalize your hearing experience with ease

Fine-tune your ReSound LiNX<sup>2</sup> hearing aids to your unique taste with the ReSound apps.

## ReSound Smart™ app:


- Adjust volume, treble and bass
- Create favorite programs and your hearing aids will switch to them in the places you visit most – or choose to have them activate automatically
- Find misplaced hearing aids with the “Find my hearing aid” functionality
- Enhance your listening comfort regardless of noise levels with the new sound enhancement features

## ReSound Relief™ app:

- The ReSound Relief app offers layered soundscapes, guided relaxation exercises, and additional activities that create a completely personalized tinnitus solution.

## ReSound Control™ app:

- Remote control access so you can adjust volume, change hearing aid programs and control your ReSound wireless accessories. All without touching – and drawing attention to – your hearing aids.



ReSound Smart app

For a complete listing of compatible phones, visit [resound.com/compatibility](https://resound.com/compatibility)



DESIGN

Fit, style and choice are all yours



New full range of hearing aid designs



Modern, discreet and durable, ReSound LiNX<sup>2</sup> hearing aids are available in a full family of new designs and colors. So when it comes to comfort and lifestyle, you're sure to find the perfect fit.

#### More comfortable

Our receiver-in-ear and behind-the-ear designs include SureFit™: a new, smaller design for the part of the hearing aid that fits into your ear canal. They may be tiny, but they are big on listening comfort and sound quality thanks to improved sound performance.

#### More durable

Our new designs stand up to the punishment of everyday life. Constructed from high-quality, durable materials, every component is coated with iSolate™ nanotech inside and out to repel moisture, dust, sweat and earwax - all substances that can cause corrosion.

# Connect, interact and engage effortlessly

Imagine stereo sound from your favorite music or movie, streamed directly from your iPhone, iPad® or iPod® touch. Make phone or video calls with your loved ones. Work out to an encouraging voice on your favorite fitness app. Enjoy it all without worrying about wires or something hanging around your neck.

ReSound LiNX<sup>2</sup> hearing aids are Made for iPhone® hearing aids, which means they connect with iPhone, iPad or iPod touch devices to make your hearing aids work like Bluetooth stereo headphones.



CONNECTIVITY

For a complete list of compatible devices go to: [www.resound.com/compatibility](http://www.resound.com/compatibility)

\* Made for iPhone connectivity is available in all behind-the-ear and receiver-in-ear styles and most custom styles.

# Connect with ReSound wireless accessories

## Stream stereo sound to your hearing aids

Don't have an iPhone? You can still stream stereo sound to your ReSound LiNX<sup>2</sup> hearing aids. Simply pair any Bluetooth-enabled phone with the ReSound Unite Phone Clip+ and use your hearing aids as headphones for phone calls, or for streaming stereo sound from videos, playlists or podcasts.

The Phone Clip+ also works as a basic remote control for your hearing aids. You can adjust the volume and mute/unmute the sound you stream to your hearing aids.



Wireless connectivity with the  
ReSound Unite Phone Clip+.



ReSound Unite TV Streamer 2

## Enjoy family movie night

With the ReSound Unite TV Streamer 2 plugged into your TV, stereo or computer, you can stream stereo sound directly into your ReSound LiNX<sup>2</sup> hearing aids. You can listen at your volume while the people around you listen at theirs.

And, because you can adjust the volume of your surroundings and the TV independently, you can still chat with family and friends while you watch TV together.



ReSound  
Micro Mic

ReSound  
Multi Mic

## Extend your listening range

The ReSound Micro Mic and ReSound Multi Mic are discreet, clip-on microphones that pair to ReSound Smart Hearing aids\* in seconds. Clip either the ReSound Micro Mic or ReSound Multi Mic on to the clothing of the person you want to talk to, and enjoy the conversation, even in background noise - and you can still hear them even if they are over 80 feet away. Plus, the ReSound Multi Mic doubles as a table microphone, so you can hear everyone sitting around it. It is also compatible with loop and FM systems\*\*, and has a line-in, so you can stream virtually any sound directly to your hearing aids.

\*Compatible with all ReSound wireless hearing aids

\*\*FM receiver required



## ReSound LiNX<sup>2</sup> delivers relief from your tinnitus

ReSound LiNX<sup>2</sup> hearing aids are equipped with a built-in Tinnitus Sound Generator (TSG). TSG provides personalized sound therapy with NatureSounds™, a new series of water-inspired sounds. And you can stream sound directly from your smartphone to your ReSound LiNX<sup>2</sup> hearing aids, which means access to any sound you find pleasing.



The new NatureSounds feature offers soothing water-inspired sounds like “ocean”, “beach surf” and “calming waves” to help provide further relief for your tinnitus patients.

When it comes to your tinnitus treatment, there’s no such thing as too much personalization. To expand the tinnitus therapy available through ReSound LiNX<sup>2</sup>, there are two apps available.

### ReSound Relief app

The ReSound Relief app offers layered soundscapes, guided relaxation exercises, and additional activities that create a completely personalized tinnitus solution.



### ReSound Smart App

The ReSound Smart app allows you to easily select and switch between the sounds that provide you the most relief from your tinnitus\*.





ReSound helps people rediscover hearing with technologies that emulate the function of the natural ear. Our aim is to give you the opportunity to live a rich, active and fulfilling life unaffected by hearing loss.

ReSound LiNX<sup>2</sup> is the new benchmark in smart hearing, giving you exceptional sound quality, speech understanding and spatial awareness. Stream stereo sound directly from your iPhone, personalize your hearing experience and engage effortlessly with family, friends and colleagues in any listening situation.

[resound.com/linx2](https://resound.com/linx2)

 [facebook.com/resoundhearing](https://facebook.com/resoundhearing)

 [twitter.com/resoundus](https://twitter.com/resoundus)

 [youtube.com/resoundus](https://youtube.com/resoundus)

### Manufacturer according to FDA:

#### ReSound North America

8001 Bloomington Freeway  
Bloomington, MN 55420  
1-888-735-4327  
[resound.com](https://resound.com)

#### ReSound Government Services

8001 Bloomington Freeway  
Bloomington, MN 55420  
1-800-392-9932  
[resound.com/veterans](https://resound.com/veterans)

### Manufacturer according to Health Canada:

#### ReSound Canada

303 Supertest Road  
Toronto, Ontario M3J 2M4  
1-888-737-6863  
[resound.com](https://resound.com)

The trademarks listed are owned and used by The GN ReSound Group and its related affiliates. © 2016. Apple, the Apple logo, iPhone, iPad and iPod touch are trademarks of Apple Inc., registered in the US and other countries. App Store is a service mark of Apple Inc. Samsung, Galaxy and Galaxy S5 are trademarks of Samsung Electronics. Android and Google Play are registered trademarks of Google Inc. The Android robot is reproduced or modified from work created and shared by Google and used according to terms described in the Creative Commons 3.0 Attribution License. Dolby, Dolby Digital and the double-D symbol are registered trademarks of Dolby Laboratories. Bluetooth is a trademark of Bluetooth SIG, Inc.

ReSound GN



MK604276

MK604276 Rev E 2016.08