

ReSound LiNX Quattro™

Confidence in the finer details



ReSound GN

Introducing
ReSound LiNX Quattro™

Our newest premium hearing aid is beyond compare with richer sound quality, longer battery life, an unrivaled full spectrum of streaming, support and personalization with ReSound Assist, and the world's most advanced rechargeable solution - all in one product.

Right down to the finer details, ReSound LiNX Quattro offers you a brilliant experience with layers of sound - and inspires confidence in the moments that matter most.



* Actual size

ReSound LiNX Quattro™



Brilliant sound experience

CLEARER SOUNDS

Built on years of innovation, ReSound LiNX Quattro utilizes the full microphone range allowing sound to be processed in a cleaner way - free of distortion.

+


FULLER SOUNDS

There are virtually no limits to your sound experience as ReSound LiNX Quattro provides you with a larger range of sounds and access to nuances in tones that were previously unavailable to you.

=

RICHER EXPERIENCE

With improved clarity and access to more sounds, coupled with a stronger signal between your hearing aids, your favorite environments become richer with ReSound LiNX Quattro.




ReSound LiNX Quattro is preferred 95% of the time for music listening, compared to other premium hearing aids.*

- ReSound LiNX Quattro
- Other premium MFi hearing aids

* Jespersen et al, 2018.

Unrivaled full spectrum of streaming

ReSound LiNX Quattro provides you with an exceptional streaming experience from almost any audio source. And with increased power, we have virtually eliminated any lag time when streaming directly into the hearing aids or through a ReSound wireless accessory.



When streaming directly from an iPhone, listeners preferred the sound quality with ReSound LiNX Quattro at least 64% of the time.*

- ReSound LiNX Quattro
- All other premium hearing aids

Find out more about the streaming experience and the full compatibility list of smartphones at resound.com





Optimize anytime, anywhere with the ReSound Smart 3D app

ReSound LiNX Quattro hearing aids automatically adjust to your environment, but you can personalize the details of your hearing experience even further. With everything you need for discreet and quick adjustments, you are in control of how you experience the sounds around you.

With ReSound Smart 3D app, you will have access to:

- Sound Enhancements
- Save Favorites and GeoTag Locations
- Find My Hearing Aids



ReSound Assist

Support and personalization

HEARING CARE AT YOUR CONVENIENCE

ReSound Assist is revolutionizing the way you receive hearing care. ReSound Assist allows you to connect directly with your hearing care professional through the ReSound Smart 3D app. ReSound Assist allows you to send requests for additional fine-tuning adjustments and get updated settings straight from your hearing care professional - without an office visit.



Find out more about ReSound Assist at
resound.com

The world's most advanced rechargeable solution

ReSound LiNX Quattro now introduces the most advanced rechargeable hearing aid – with a lithium-ion battery that is completely sealed for protection. With its small, discreet size and long-lasting battery, you get power that lasts more than a full day – without having to fiddle with any batteries.


Battery status is easy to monitor with the ReSound Smart 3D app and the LED lights on the charging case indicate the power status for both your hearing aids and the charger.

FULL CONFIDENCE WITH A COMPLETE CHARGE

At night, simply place your hearing aids in the inductive, easy-to-use charging case to prepare for the next day. If you do forget to charge, the case's built in battery holds three full charges of power, giving you the freedom to charge anywhere.

CHARGING TIME	HOURS OF BATTERY LIFE
3 hours =	30 24 with 50% streaming time
1 hour =	16
30 minutes =	8

Find out more about rechargeability at [resound.com](https://www.resound.com)



High performance style options

The same clearer, fuller and richer sound is available in a variety of hearing aid styles. In addition to the small and long-lasting rechargeable solution, ReSound LiNX Quattro offers a replaceable battery solution and a family of all wireless custom models. With low power consumption and unparalleled sound quality, you can enjoy the serenity of nature or bustle of the city with complete confidence.

CUSTOM STYLES



CUSTOM COLORS



RIE COLORS



RIE STYLES



WORLD'S MOST ADVANCED RECHARGEABLE SOLUTION



REPLACEABLE BATTERY SOLUTION

People with hearing loss are at the heart of everything we do because we know they're placing their trust in us. This is why we've been pioneering innovative new technologies, with consistent quality, for over 150 years to help people feel more involved, connected and in control, transforming their lives through the power of sound. Available exclusively through audiologists and hearing care professionals.

ReSound is part of the GN Group – pioneering great sound from world-leading ReSound hearing aids to Jabra office headsets and sports headphones. Founded in 1869, employing over 5,000 people, and listed on NASDAQ OMX Copenhagen, GN makes life sound better.

Find out more about how ReSound LiNX Quattro™ can provide confidence in the finer details at resound.com.

 facebook.com/resoundhearing

 twitter.com/resoundus

 youtube.com/resoundus

Manufacturer according to FDA:

GN ReSound North America

8001 E Bloomington Freeway
Bloomington, MN 55420
USA
1-888-735-4327
resound.com

Manufacturer according to Health Canada:

ReSound Canada

303 Supertest Road
Toronto, ON M3J 2M4
Canada
1-888-737-6863
resound.com

ReSound Government Services

8001 E Bloomington Freeway
Bloomington, MN 55420
USA
1-800-392-9932
resound.com/veterans

© 2019 GN Hearing A/S. All rights reserved. ReSound is a trademark of GN Hearing A/S. Apple, the Apple logo, iPhone, iPad and iPod touch are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc., registered in the U.S. and other countries. Google Play and the Google Play logo are trademarks of Google Inc.



MK605091

MK605091 Rev B 2019.07

ReSound GN
.....