

BILL TO

Acc't #:

Office:

Address:

Address:

City/State:

Zip: P.O.:

SHIP TO

Acc't #:

Date: Phone:

Contact name:

Email:

Facility:

Address:

Address:

City/State:

Zip:

1

PATIENT DATA

Patient's name:

SSN:

LAST:

FIRST:

Previous user: ☐ YES ☐ NO

Audiogram data:

250 500 750 1k 1.5k 2k 3k 4k 6k 8k

Left air:

Bone:

LDL:

Right air:

Bone:

LDL:

2

SPECIAL INSTRUCTIONS

(PLEASE PRINT CLEARLY)

3

MODEL AND OPTIONS

2.4GHz wireless
ReSound LiNX Quattro
Portable charging case included in order of RE961

2.4GHz wireless
ReSound LiNX 3D

Mini RIC	Mini RIC	Mini RIC	RIC
L R	L R	L R	L R
RE961	—	—	RE962
☐ ☐			☐ ☐
—	LT961	LT961-Z	LT962
	☐ ☐	☐ ☐	☐ ☐
Battery	Lithium-ion rechargeable	#312	#312
Program selector	Programmable Push button	Programmable Push button	Multifunction Switch
Volume control	Programmable Push button	Programmable Push button	Multifunction Switch
RIC receivers (power level)			
Low power (LP)	☐ ☐	☐ ☐	☐ ☐
Medium power (MP)	☐ ☐	☐ ☐	☐ ☐
High power (HP)	☐ ☐	☐ ☐	☐ ☐
Ultra power (UP)*	☐ ☐	☐ ☐	☐ ☐
Color (keychain color)			
Dark brown (4)	☐ ☐	☐ ☐	☐ ☐
Beige (5)	☐ ☐	☐ ☐	☐ ☐
Light blonde (6)	☐ ☐	☐ ☐	☐ ☐
Medium blonde (7)	☐ ☐	☐ ☐	☐ ☐
Monza red (8)	☐ ☐	☐ ☐	☐ ☐
Black (9)	☐ ☐	☐ ☐	☐ ☐
Anthracite (10)	☐ ☐	☐ ☐	☐ ☐
Silver (11)	☐ ☐	☐ ☐	☐ ☐
Forest camo	☐ ☐	☐ ☐	☐ ☐
Desert camo	☐ ☐	☐ ☐	☐ ☐
Ocean camo	☐ ☐	☐ ☐	☐ ☐
Pearl white (12)	☐ ☐	☐ ☐	☐ ☐
Ocean blue (13)	☐ ☐	☐ ☐	☐ ☐
Sterling gray (14)	☐ ☐	☐ ☐	☐ ☐
Gloss black (15)	☐ ☐	☐ ☐	☐ ☐
Gloss anthracite (16)	☐ ☐	☐ ☐	☐ ☐
Gloss medium blonde (17)	☐ ☐	☐ ☐	☐ ☐
White	☐ ☐	☐ ☐	☐ ☐

* UP receiver only available in Encased earmolds
(#) Color keychain reference

☐ AVAILABLE ☐ DEFAULT

4a

RESOUND RIC RECEIVERS AND DOMES

RIC Receiver tube length		DOMES		Size	STD	POWER
L	R	L	R			
0	☐ ☐	Large	☐ ☐	Large	☐ ☐	
1	☐ ☐	Medium	☐ ☐	Medium	☐ ☐	
2	☐ ☐	Small	☐ ☐	Small	☐ ☐	
3	☐ ☐	Tulip	☐ ☐			
4	☐ ☐			Large	☐ ☐	
				Medium	☐ ☐	
				Small	☐ ☐	

OR

4b

CUSTOM RIC EARMOLD OPTIONS



For proper fit, measure from the top of the ear to the **center** of the ear canal using the ReSound measuring tool.

- ☐ (\$) ReSound Remote Control 2
- ☐ (\$) ReSound TV Streamer 2
- ☐ (\$) ReSound Phone Clip+
- ☐ (\$) ReSound Multi Mic
- ☐ (\$) ReSound Micro Mic
- ☐ (\$) ReSound LiNX Quattro RE961
Portable Charging Case (additional)

(\$)
Additional
charge for
wireless
accessories

MATERIAL

- ☐ Hard (acrylic) ☐ Soft (silicone)
(n/a for Encased and Hollow Cavity)

ENCASED



MICROMOLD



HOLLOW CAVITY



SKELETON




INSTRUMENT INFORMATION

MODEL _____

TRUFIT™ IMPRESSION—THE 16/4 RULE

Take an **OPEN JAW** impression when:

- Ear geometry lacks retention
- Patient has severe TMJ movement
- Instrument migrates out of ear
- Instrument is loose or has feedback



RECEIVER		Include (check): YES NO		Size: _____	
Low power (LP)		L	R		
Medium power (MP)		L	R		
High power (HP)		L	R		
Ultra power (UP) (Encased only)		L	R		
SHELL COLOR					
Clear		L	R		
Light		L	R		
Medium		L	R		
Dark		L	R		
Rose (n/a for Encased nor Hollow Cavity)		L	R		
EarLusion Light (n/a for Encased nor Hollow Cavity)		L	R		
Espresso (hard only)		L	R		
Red/Blue		L	R		
FACEPLATE COLOR (Encased only)					
Light		L	R		
Beige		L	R		
Medium		L	R		
Dark		L	R		
Espresso		L	R		
Anthracite		L	R		
Clear		L	R		
CANAL LENGTH					
Factory select		L	R		
As marked		L	R		
VENTING					
Factory select		L	R		
MOV (Semi-IROS vent modification recommended)		L	R		
SAV		L	R		
Pressure		L	R		
None		L	R		
VENT MODIFICATION					
Semi-IROS		L	R		
IROS (n/a for Hollow Cavity)		L	R		
WAX PROTECTION (Encased and hard only, n/a for Hollow Cavity)					
HF3		L	R		
CeruSTOP (default for Encased)		L	R		
None (default for hard, STD for Hollow Cavity)		L	R		
OTHER OPTIONS					
Removal cord		L	R		
Blue/Red dots		Size (check): SMALL or LARGE			
Patient initials	_____				
RETENTION					
Canal Lock (n/a for Skeleton)		L	R		
Helix Lock (n/a for Hollow Cavity and Skeleton)		L	R		
Skeleton Lock (n/a for Skeleton)		L	R		
Semi-Skeleton Lock (n/a on Hollow Cavity and Skeleton)		L	R		
PLEASE SEND ☐ Air bills ☐ Impression mailers					

AVAILABLE DEFAULT

